

RF hand pieces for Face Skin Tightening and Body Contouring

Dr. Betty Czajkowsky, M.D.

SharpLight Technologies Medical Advisor

INTRODUCTION

Radiofrequency dermal heating devices have recently emerged as an effective, non-invasive, aesthetic treatment modality to improve the appearance of cellulite, to eliminate face wrinkles and rhytids, to tighten lax skin and reduce abdominal and thighs circumference. Radiofrequency devices heat under the skin by passing rapidly oscillating electrical currents through the tissues. In monopolar devices current is passed between a small electrode in contact with the treatment area and a large return electrode pad, usually placed under the buttocks, while with bipolar devices the current is passed between two electrodes in contact with the treatment area. Bipolar devices offer better control over the treated tissue volume and require lower powers to heat up a specific volume of under skin. With bipolar devices skin heating depth is a function of the distance between the electrodes, the RF power and the dwell time over the same tissue area. By proper adjustment of these parameters both the dermis and subcutaneous layers can be heated to therapeutic levels while preserving the epidermis. Treatments can be performed with absolutely no downtime and can often result in dramatic changes in body contours as well as face skin texture. Face, neck, arms, abdomen, buttocks and thighs are all amenable to RF treatments.

Radiofrequency aesthetic face skin tightening treatments have been documented to be extremely safe with only minimal and transient side effects such as erythema and slight edema. To maintain this level of safety it was found that anesthesia should be avoided to allow patient feedback regarding pain and heat sensation, which is totally bearable and even pleasant.

Heat produced by RF in the dermis and subcutaneous layer causes immediate tightening of collagen and induces new collagen production. Increase in local blood circulation, promotes fibrous tissue breakdown, drainage of fatty deposits through the lymphatic system and disintegration of fat/cellulite cells. Adipose tissues “mobilize” or release stored fat to make it available as energy for metabolic activities by a physiological mechanism called lipolysis. Like any chemical reaction all enzyme-mediated biological reaction rates, including those of lipolysis, are temperature sensitive. Increasing the temperature of the fat cells by 10°C increases the lipolytic process rate by factor of 2. Increasing the fat cells temperature, above 43°C, results in protein denaturation and irreversible damage. All these mechanisms contribute to the skin tightening effects, skin texture improvements, and reduction in the appearance of cellulite and circumference reduction, all of which have been demonstrated in clinical studies with various RF devices. In these studies facial areas usually required 4-8 weekly sessions while body areas require a course of 8-12 weekly sessions.

Patients usually feel the skin tightening effect immediately following each treatment however the full result becomes evident following the end of the treatment sessions due to the build up of new collagen in the area. Results are expected to last 1-2 years and maintenance treatments can be continued on a regular basis or repeated as desired by the patient.

Finally, RF aesthetic treatments can be combined with a variety of body contouring and skin texture complimentary techniques ranging from cosmetology external products, tissue massage techniques such as LPG, different dietary and physical training programs and other phototherapy technologies as Infrared and Lasers and even as a maintenance tool after plastic surgery.

MATERIALS AND METHODS

A new bipolar RF hand pieces for the Formax Plus and Omnimax Aesthetic Stations has been tested on twenty (20) volunteers, 14 females, 6 males, ages between 25 and 45, with face lax skin and various degrees of fat/cellulite in the abdomen and thighs and stretch marks in the buttocks. The twenty (20) volunteers received eight (8) treatments, once a week. This handpiece emits continuous RF powers of up to 50 Watts. Time of the treatment sessions vary depending on the area treated: face and neck around 20 minutes and body areas like abdomen and both legs from 40 to 60

minutes. A non-contact infrared thermometer was used to monitor external skin therapeutic temperature during sessions in order to mantain a constant over 40°C during all the treatment session. Standardized photographs were taken at baseline, before each weekly treatment and at the 1 month follow-up visit. Patient satisfaction scores were also recorded at each visit. Finally circumference measurements were taken on abdomen and thighs at baseline and at the final follow-up visit. Each session's data was file for posterior comparative objective evidence.



Redness and heat sensation lasting for about 30 minutes after the session

Treatment contraindications include:

Patients with any active "metal" implants, such as a cardiac pacemaker, treatment over permanent implants, current or history of cancer, diabetes, congestive heart disease, epilepsy, neurological disorders or active infections, significant skin irregularities in treatment area, including dry skin, open wounds, major scars, microvascular disorders, or any inflammatory skin conditions, any recent intervention in the target area which requires a waiting period, use of oral isotretinoin within two

months, or use of topical isotretinoin within one month of treatment, immunosuppressive diseases, including AIDS and HIV infection, or use of immunosuppressive medications, history of bleeding coagulopathies or use of anticoagulants, history of keloid scarring, skin atrophies, or poor wound healing, adverse reaction to the RF conductive lotion, pregnancy and/ or nursing

The circular, bi-polar hand piece used in this study operates at a frequency of 1 MHz and a diameter of 7.5 cm.



FACE hand piece



BODY hand piece

RESULTS

At 1 month after the last treatment, 90% of the subjects exhibited moderate to significant improvement in skin laxity and skin texture on the face, neck, arms, abdomen, buttocks and thighs. A significant reduction in abdominal and thigh circumference was found following RF treatment.

Areas with cellulite showed marked improvement. Patients reported high

overall satisfaction with the treatment. No serious complications were recorded.

- Average abdominal circumference reduction in 8 treatments was 4 to 5 cm
- Average tight (back legs) reduction in 8 treatments was 3 to 4 cm
- Face and neck had a positive subjective evaluation after 8 treatments

Patient satisfaction scores:

- Very satisfied 80 %

- Satisfied 15 %

- Not satisfied: 5 %

CONCLUSION

The new bipolar RF hand pieces for Sharplight's Aesthetic Stations has been demonstrated to be highly effective and safe for non-invasive skin tightening and body shaping. Quantitative and qualitative results were documented in the treatment of face and neck skin laxity and of cellulite and circumference reduction of abdomen, thighs and buttocks. Very high patient satisfaction was achieved. The procedure is easy to apply, involves no discomfort to the patient and requires no downtime whatsoever.

Clearly visible results were obtained after eight weekly sessions only, with improvement continuing after treatments are terminated. Further treatment sessions or periodic maintenance sessions may be administered as required or desired by each individual patient.

The availability of this new RF hand pieces significantly enhances the capabilities of the SLT Aesthetic Stations, broadening its clinical indications to include skin tightening on all body areas.

B/A PHOTO

Before



After 4 sessions



Skin tightening and neocollagenesis

Before



After 8 sessions



Lost of 6 cms circumference reduction

Before



After 8 sessions



Lost of 4 cms circumference reduction

Fat reduction, cellulite reduction,

Stretch marks reduction

REFERENCES

Clinical and histopathological results following TriPollar radiofrequency skin treatments. **Kaplan H, Gat A.** *J Cosmet Laser Ther.* 2009 Jun;11(2):78-84.

Nonsurgical nonablative treatment of aging skin: radiofrequency technologies between aggressive marketing and evidence-based efficacy.

Atiyeh BS, Dibo SA. *Aesthetic Plast Surg.* 2009 May;33(3):283-94. Epub 2009 May 13.

Monopolar radiofrequency treatment of human eyelids: a prospective, multicenter, efficacy trial. **Biesman BS, Baker SS, Carruthers J, Silva HL, Holloman EL.** *Lasers Surg Med.* 2006 Dec;38(10):890-8.

Is topical anesthesia useful in noninvasive skin tightening using radiofrequency?

Kushikata N, Negishi K, Tezuka Y, Takeuchi K, Wakamatsu S. *Dermatol Surg.* 2005 May;31(5):526-33.

Non-ablative skin tightening with radiofrequency in Asian skin.

Kushikata N, Negishi K, Tezuka Y, Takeuchi K, Wakamatsu S. *Lasers Surg Med.* 2005 Feb;36(2):92-7.

Selective electro-thermolysis in aesthetic medicine: a review.

Sadick NS, Makino Y. *Lasers Surg Med.* 2004;34(2):91-7.

Clinical applications of radiofrequency: nonsurgical skin tightening

Hodgkinson DJ. *Clin Plast Surg.* 2009 Apr;36(2):261-8, viii.

The nonablative radiofrequency for rejuvenation. **Sukal SA, Geronemus RG.** *Clin Dermatol.* 2008 Nov-Dec;26(6):602-7.

Tissue tightening: a hot topic utilizing deep dermal heating. **Gold MH.** *J Drugs Dermatol.* 2007 Dec;6(12):1238-42

Clinical, laboratory, and MRI analysis of cellulite treatment with a unipolar radiofrequency device. **Goldberg DJ, Fazeli A, Berlin AL.**

Dermatol Surg. 2008 Feb;34(2):204-9; discussion 209. Epub 2007 Dec 17.

Bipolar radiofrequency in the treatment of dermatologic imperfections: clinicopathological and immunohistochemical aspects.

Montesi G, Calvieri S, Balzani A, Gold MH. *J Drugs Dermatol.* 2007 Sep;6(9):890-6.

Nonablative cutaneous remodeling using radiofrequency devices.

Alster TS, Lupton JR. *Clin Dermatol.* 2007 Sep-Oct;25(5):487-91.

Nonsurgical tightening of skin laxity: a new radiofrequency approach.

Rusciani A, Curinga G, Menichini G, Alfano C, Rusciani L. *J Drugs Dermatol.*

2007 Apr;6(4):381-6.

Monopolar radiofrequency skin tightening.

Abraham MT, Mashkevich G. *Facial Plast Surg Clin North Am.* 2007 May;15(2):169-77, vi.

Skin tightening with a combined unipolar and bipolar radiofrequency device.

Mayoral FA. *J Drugs Dermatol.* 2007 Feb;6(2):212-5.

The role of deep heating for noninvasive skin rejuvenation

Dierickx CC. *Lasers Surg Med.* 2006 Oct;38(9):799-807.

Monopolar radiofrequency facial tightening: a retrospective analysis of efficacy and safety in over 600 treatments.

Weiss RA, Weiss MA, Munavalli G, Beasley KL. *J Drugs Dermatol.* 2006 Sep;5(8):707-12.

Monopolar radiofrequency.

Burns JA. *Aesthet Surg J.* 2005 Nov-Dec;25(6):638-42.

Nonablative radiofrequency treatment of facial laxity.

Fisher GH, Jacobson LG, Bernstein LJ, Kim KH, Geronemus RG. *Dermatol Surg.* 2005 Sep;31(9 Pt 2):1237-41; discussion 1241.

Multipass vector (mpave) technique with nonablative radiofrequency to treat facial and neck laxity. **Finzi E, Spangler A.** *Dermatol Surg.* 2005 Aug;31(8 Pt 1):916-22.

Current concepts in nonablative radiofrequency rejuvenation of the lower face and neck. **Abraham MT, Vic Ross E.** *Facial Plast Surg.* 2005 Feb;21(1):65-73.

Radiofrequency nonablative tissue tightening. **Koch RJ.**

Facial Plast Surg Clin North Am. 2004 Aug;12(3):339-46, vi.

Improvement of neck and cheek laxity with a nonablative radiofrequency device: a lifting experience.

Alster TS, Tanzi E. *Dermatol Surg.* 2004 Nov;30(11):1430; author reply 1430.

Non-surgical radiofrequency facelift.

Narins DJ, Narins RS. *J Drugs Dermatol.* 2003 Oct;2(5):495-500.

The medical face lift: a noninvasive, nonsurgical approach to tissue tightening in facial skin using nonablative radiofrequency.

Ruiz-Esparza J, Gomez JB. *Dermatol Surg.* 2003 Apr;29(4):325-32; discussion 332.

Circumference reduction and cellulite treatment with a TriPollar radiofrequency device: a pilot study.

Manuskiatti W, Wachirakaphan C, Lektrakul N, Varothai S. J Eur Acad Dermatol Venereol. 2009 Jul;23(7):820-7. Epub 2009 Apr 8.

Histological findings in adipocytes when cellulite is treated with a variable-emission radiofrequency system.

Trelles MA, van der Lugt C, Mordon S, Ribé A, Al-Zarouni M. Lasers Med Sci. 2009 Mar 26. [Epub ahead of print]

Adipocyte membrane lysis observed after cellulite treatment is performed with radiofrequency. **Trelles MA, Mordon SR. : Aesthetic Plast Surg.** 2009 Jan;33(1):125-8. Epub 2008 Nov 27.

Unipolar radiofrequency treatment to improve the appearance of cellulite.

Alexiades-Armenakas M, Dover JS, Arndt KA. J Cosmet Laser Ther. 2008 Sep;10(3):148-53.

An evidence-based assessment of treatments for cellulite.

Wanner M, Avram M. J Drugs Dermatol. 2008 Apr;7(4):341-5.

Effect of controlled volumetric tissue heating with radiofrequency on cellulite and the subcutaneous tissue of the buttocks and thighs.

Emilia del Pino M, Rosado RH, Azuela A, Graciela Guzmán M, Argüelles D, Rodríguez C, Rosado GM. J Drugs Dermatol. 2006 Sep;5(8):714-22.

SharpLight™

SharpLight Technologies Ltd. 33
Lazarov St. Rishon Le Zion, Israel
75654
Tel: + 972 (0) 3 961 1969
Fax: + 972 (0) 3 952 8890
info@sharplighttech.com
www.sharplighttech.com
MRK-00026 Rev 2